

New Post-Evisceration Steep Radius Conformers


Products of Quality Without Compromise



Now you can finally match the conformer radius to the globe shape when performing an evisceration procedure. No more need to worry about central blanching that may cause your closure to become exposed possibly leading to an extrusion of the implant.

Look at this side view of a steep radius evisceration conformer resting on an ocular prosthesis with a radius of 44.50 diopters.

Notice how nicely the conformer fits the shape of the eye, which is very similar to that of an eviscerated globe.



Our new steep radius evisceration conformers are just the ticket when you have concerns about possible pressure on the surgical closure following an evisceration or in some cases an enucleation.

We offer three sizes, all packaged in sterile pouches. The sizes are: Small: SRC-01, 21x18mm, Medium: SRC-02, 22x21mm and Large: SRC-03, 24x22mm. You can purchase them one at a time or by the box. Order a box of 10 all the same size, or choose our assorted box of 10: SRC-10-PA with 2 small, 4 medium and 4 large.

Kolberg Ocular Supplies, Inc.

Call toll free: 866-259-5004 ~ Fax: 866-950-4704

Email: corneashields@aol.com

Web: www.corneaprotectors.com